# Package 'sleepr'

January 9, 2019

January 9, 2019
Title Analyse Activity and Sleep Behaviour
<b>Date</b> 2018-10-04
Version 0.3.0
Description Use behavioural variables to score activity and infer sleep from bouts of immobility. It is primarily designed to score sleep in fruit flies from Drosophila Activity Monitor (TriKinetics) and Ethoscope data.  It implements sleep scoring using the ``five-minute rule" (Hendricks et al. (2000) <doi:10.1016 s0896-6273(00)80877-6="">), activity classification for Ethoscopes (Geissmann et al. (2017) <doi:10.1371 journal.pbio.2003026="">) and a new algorithm to detect when animals are dead.</doi:10.1371></doi:10.1016>
<b>Depends</b> R (>= 3.00), behave
Imports data.table
Suggests testthat, covr, knitr
License GPL-3
Encoding UTF-8
LazyData true
<pre>URL https://github.com/rethomics/sleepr</pre>
BugReports https://github.com/rethomics/sleepr/issues RoxygenNote 6.1.0 Roxygen list(markdown = TRUE)
R topics documented:
bout_analysis

bout\_analysis

Index 8

bout\_analysis

Find "bouts" in categorical time series

# **Description**

This function is used to find contiguous regions of unique value in a – potentially irregular/heterogeneous – univariate categorical time series.

# Usage

```
bout_analysis(var, data)
```

# Arguments

var name of the variable to use from data

data data.table containing behavioural variable from or one multiple animals. When

it has a key, unique values, are assumed to represent unique individuals (e.g. in a behavr table). Otherwise, it analysis the data as coming from a single animal.

data must have a column t representing time.

#### Value

an object of the same type as data (i.e. data.table::data.table or behavr::behavr). Each row is a specific bout characterised by three columns.

- t its *onset*
- duration its length
- <var> a column with the same name as var. The value of var for this bout.

### References

• The relevant rethomic tutorial section – on sleep analysis

# See Also

- sleep\_annotation to generate a binary sleep variable
- rle run length encoding function on which this analysis is based

curate\_dead\_animals 3

#### **Examples**

curate\_dead\_animals

Remove - irrelevant - data after individual have died

# **Description**

This function detects when individuals have died based on their first (very) long bout of immobility. Following data (which may include spurious artefact of movement) are removed.

#### Usage

```
curate_dead_animals(data, moving_var = moving, time_window = hours(24),
   prop_immobile = 0.01, resolution = 24)
```

#### **Arguments**

data.table containing behavioural variable from or one multiple animals. When it has a key, unique values, are assumed to represent unique individuals (e.g. in a behavr table). Otherwise, it analysis the data as coming from a single animal. data must have a column t representing time.
logical variable in data used to define the moving (alive) state (default is moving)
window during which to define death (default is one day)
proportion of immobility that counts as "dead" during time_window (see details)
how much scanning windows overlap. Expressed as a factor (see details).

#### **Details**

This function scans the time series looking for "death" in the right (future) data, within time\_window. Death is defined as mean(moving\_var) < prop\_immobile within a time window. Moving window start every time\_window/resolution. resolution = 1 is fast but means no overlap. The default would score an animal as dead it does not move more than *one percent of the time* for at least *one day*. All data following a "death" event are removed.

4 motion\_detectors

#### Value

an object of the same type as data (i.e. data.table::data.table or behavr::behavr).

# See Also

• sleep\_annotation – to score movement and slepe

# **Examples**

```
dt1 <- toy_activity_data()</pre>
#all movement after day 3 is set at 0
dt1[t > days(3), moving := FALSE]
# one artefact of movement is detected at day 3.5
dt1[t == days(3.5), moving := TRUE]
dt2 <- curate_dead_animals(dt1)</pre>
dt3 <- curate_dead_animals(dt1,prop_immobile = 0)</pre>
## Not run:
library(ggplot2)
ggplot(data=dt1[,test:=1],aes(t, as.numeric(moving))) +
  geom_line(data=dt1[,test:=1]) +
  geom_line(data=dt2[, test:=2])+
  geom_line(data=dt3[, test:=3])+
  facet_grid(test ~ .)+
  scale_x_time()
## End(Not run)
```

motion\_detectors

Motion detector for Ethocope data

#### **Description**

Defines whether a single animal is moving according to:

## Usage

```
max_velocity_detector(data, time_window_length,
    velocity_correction_coef = 0.003, masking_duration = 6)
max_velocity_detector_legacy(data, velocity_threshold = 0.006)
virtual_beam_cross_detector(data, time_window_length)
```

motion\_detectors 5

#### **Arguments**

data

data.table::data.table containing behavioural variables of *a single animal* (no id). It must have the columns xy\_dist\_log10x1000(for computing subpixel velocity), x(beam cross), t and has\_interacted (whether a stimulus was delivered).

time\_window\_length

number of seconds to be used by the motion classifier. This corresponds to the sampling period of the output data.

velocity\_correction\_coef

an empirical coefficient to correct velocity with respect to variable framerate.

masking\_duration

number of seconds during which any movement is ignored (velocity is set to 0) after a stimulus is delivered (a.k.a. interaction).

velocity\_threshold

uncorrected velocity above which an animal is classified as 'moving' (for the legacy version).

#### **Details**

- Validated and corrected subpixel velocity (max\_velocity\_detector), the most rigorous
- Uncorrected subpixel velocity (max\_velocity\_detector\_legacy)
- Crossing a virtual beam in the middle of the region of interest (virtual\_beam\_cross\_detector)

max\_velocity\_detector is the default movement classification for real-time ethoscope experiments. It is benchmarked against human-generated ground truth.

These functions are rarely called directly, but typically used is in the context of sleep\_annotation.

# Value

an object of the same type as data (i.e. data.table::data.table or behavr::behavr) with additional columns:

- moving Logical, TRUE iff. motion was detected.
- beam\_crosses The number of beam crosses (when the animal crosses x = 0.5 that is the midpoint of the region of interest) within the time window
- max\_velocity The maximal velocity within the time window. The resulting data is sampled at a period equals to time\_window\_length.

#### See Also

• sleep\_annotation – which requieres a motion detector

6 sleep\_annotation

sleep\_annotation

Score sleep behaviour from immobility

### **Description**

This function first uses a motion classifier to decide whether an animal is moving during a given time window. Then, it defines sleep as contiguous immobility for a minimum duration.

#### **Usage**

```
sleep_annotation(data, time_window_length = 10,
    min_time_immobile = 300, motion_detector_FUN = max_velocity_detector,
    ...)
sleep_dam_annotation(data, min_time_immobile = 300)
```

### Arguments

data

data.table containing behavioural variable from or one multiple animals. When it has a key, unique values, are assumed to represent unique individuals (e.g. in a behavr table). Otherwise, it analysis the data as coming from a single animal. data must have a column t representing time.

time\_window\_length

number of seconds to be used by the motion classifier. This corresponds to the sampling period of the output data.

min\_time\_immobile

Minimal duration (in s) of a sleep bout. Immobility bouts longer or equal to this value are considered as sleep.

motion\_detector\_FUN

function used to classify movement

... extra arguments to be passed to motion\_classifier\_FUN.

#### **Details**

The default time\_window\_length is 300 seconds – it is also known as the "5-minute rule". sleep\_annotation is typically used for ethoscope data, whilst sleep\_dam\_annotation only works on DAM2 data. These functions are *rarely used directly*, but rather passed as an argument to a data loading function, so that analysis can be performed on the go.

#### Value

a behave table similar to data with additional variables/annotations (i.e. moving and asleep). The resulting data will only have one data point every time\_window\_length seconds.

#### References

• The relevant rethomic tutorial section – on sleep analysis

sleep\_annotation 7

#### See Also

- motion\_detectors options for the motion\_detector\_FUN argument
- bout\_analysis to further analyse sleep bouts in terms of onset and length

#### **Examples**

```
dt_one_animal <- toy_ethoscope_data(seed=2)</pre>
###### Ethoscope, corrected velocity classification ########
sleep_dt <- sleep_annotation(dt_one_animal, masking_duration=0)</pre>
print(sleep_dt)
# We could make a sleep `barecode'
## Not run:
library(ggplot2)
ggplot(sleep_dt, aes(t,y="Animal 1",fill=asleep)) +
                                    geom_tile() + scale_x_time()
## End(Not run)
###### Ethoscope, virutal beam cross classification ########
sleep_dt2 <- sleep_annotation(dt_one_animal,</pre>
                             motion_detector_FUN=virtual_beam_cross_detector)
## Not run:
library(ggplot2)
ggplot(sleep_dt2, aes(t,y="Animal 1",fill=asleep)) +
                                    geom_tile() + scale_x_time()
## End(Not run)
###### DAM data, de facto beam cross classification ######
dt_one_animal <- toy_dam_data(seed=7)</pre>
sleep_dt <- sleep_dam_annotation(dt_one_animal)</pre>
## Not run:
library(ggplot2)
ggplot(sleep_dt, aes(t,y="Animal 1",fill=asleep)) +
                                    geom_tile() + scale_x_time()
## End(Not run)
```

# **Index**

```
behavr, 2, 3, 6
behavr::behavr, 2, 4, 5
bout_analysis, 2, 7
curate\_dead\_animals, 3
data.table, 2, 3, 6
data.table::data.table, 2, 4, 5
{\tt max\_velocity\_detector}, {\tt 5}
max_velocity_detector
         (motion_detectors), 4
max_velocity_detector_legacy, 5
{\tt max\_velocity\_detector\_legacy}
         (motion_detectors), 4
motion_detectors, 4, 7
rle, 2
sleep_annotation, 2, 4, 5, 6
sleep\_dam\_annotation
        (sleep_annotation), 6
virtual_beam_cross_detector, 5
virtual_beam_cross_detector
        (motion_detectors), 4
```